


Contact Quality



Sampling Rate: 128  
Battery: [Empty]

Marker

Serial Port: COM1 [Refresh]  
[Configure Marker ...]

Log Message:  
Marker(2) is inserted  
Marker(1) is inserted  
Marker(2) is inserted  
Marker(1) is inserted  
Marker(2) is inserted  
Marker(1) is inserted  
Marker(2) is inserted  
[Clear]

[Save Data]  
[Load Data]

EEG FFT Gyro Data Packets

Channel: O2  
amax: 80  
amin: -60  
fmax: 64  
fmin: 0  
length: 512  
step: 1  
dB: 10  
Windowing: HANNING

pmax: 2.50  
pmin: 0.00  
[Auto Scale]  
Custom Band  
lower: 30 Hz  
upper: 45 Hz

